## Ultimate Recover Back

cold compression therapy  
for sports injuries

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### Contains

Nylon, Polyester, Polycarbonate, Polypropylene,  
Polyethylene, Polyvinyl chloride  
**Gelpack**  
Propylene glycol, Sodium carboxymethyl cellulose

Read the instructions for use before using the product.  
Lees voor gebruik van dit medisch hulpmiddel de  
gebruiksaanwijzing. Veuillez lire attentivement la notice  
d'utilisation avant d'utiliser ce produit. Lesen Sie die  
Gebrauchsanweisung sorgfältig durch bevor Sie dieses  
Produkt verwenden.

**GB** Incorrect application can cause injuries. Consult a  
qualified person for the correct application. Stop using  
the product immediately when skin irritation occurs.  
Keep out of reach of children. Keep away from heat,  
hot surfaces, sparks, open flames and other ignition  
sources. No smoking. Store in a well ventilated place.  
Dispose of contents / container in accordance with  
local / regional / national / international regulations.

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
### BACK

#### PARTS

- 1 compressor-cover
- 2 gelpack
- 3 one strap
- 4 pump with tube
- 5 vent button

#### INSTRUCTION FOR USE

1. Put the gelpack (2) in the fridge / freezer at least 2 hours before use.
2. Put the gelpack (2) with its back (velcro tape) on the compressor-cover (1).
3. Pull one end of the strap (3) with the black Velcro, through the square hook (from IN to OUT) and fasten (image 1).
4. Pull the other end of the one strap (3) through the other end of the square hook.
5. Then adjust both sides of the Velcro tape until comfortable (image 2 & 3).
6. Connect the tube (4) to the ON/OFF switch, screw it on the two-way valve.
7. Turn the switch in the ON position.
8. Pump it up (for compression) and increase the pressure around the affected back\*.
9. Turn the switch to the OFF position and just twist and unwind the tube.
10. Brace on for 5-10 minutes, Brace off for 5-10 minutes and repeat.
11. Empty the compressor-cover before taking the brace off: Push in the vent button (5).

\*Avoid over-compression; excessive pressure may cause blood circulation disorder.

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