GB INSTRUCTION FOR USE

Read the instruction for use before using the product. Keep out of reach of children. Stop using the product immediately when skin irritation occurs. Test pad is recommended for those with high allergic tendencies. Use only on intact skin. For single use only. Store at room temperature. Avoid storing product in direct sunlight and/or high temperatures.

1. CONTENTS

1x Instruction for use
1x Kinesio tape

CONTAINS

Acrylic-polymer, cotton and polyurethane.

2. DESCRIPTION

Athlete Kinesio Tape is an elastic tape for prevention and treatment of sport injuries and sport-related aches and pains. Kinesio Tape can be used for support, pain relief and recovery of muscles and joints during or after training.

FEATURES

- effective on its own, but can be combined with other therapeutic modalities for even greater impact
- allows athletes to remain active while injured
- skin ventilation is possible through high porosity

ADVICE

Use U-Sport Pre Tape and/or U-Sport Underwrap for better bonding/adhesion. Remove the tape easily and effectively with U-Sport Tape Remover.

3. INTENDED USE / INDICATIONS FOR USE

Athlete Kinesio Tape provides excellent comfort for muscles and joints. The product is intended for use:

- only on intact skin
- people over 12 years of age
- for single use only

4. INSTRUCTIONS FOR USE

The incorrect application of kinesio tape can cause injuries. Consult a qualified person for the correct application. Storage temperatures can have influence on the adhesive strength of kinesio-tape. Can be used for a wide variety of areas of the body. Usually lasts max. 3 days per application.

INSTRUCTIONS

- Thorough cleaning and drying of the skin is essential for the proper adhesion of all (sports) tapes.
- check if the skin is clean, dry and free from any creams or lotions
 - if body hair is limiting adhesion you may need to trim/shave the area
 - elongate the muscles or joints, prior application of the adhesive bandage
 - apply the adhesive bandage with little or no stretch
 - apply with correct application method (I-shape, Y-shape, X-shape and Fork shape), see image
 - wait at least 1 hour after application before engaging an activity that causes sweating
 - do not swim or shower for at least 1 hour after application
 - dry the tape pat gently with a towel

DO NOT


Apply it very tightly because of the risk of cutting off blood circulation.

REMOVAL

- avoid lifting the skin, never rip or pull off tape forcefully
- pull with the direction of hair growth and
- while removing the tape: immediately move the thumb (with a little pressure) to keep the skin down

CAUTION

Test pad is recommended for those with allergic tendencies. Stop using the product immediately when skin irritation occurs

5. CONTRAINDICATIONS

DO NOT USE

- if sensitive to one of the ingredients or skin irritation occurs
- history of hyper-sensitive skin or pre-existing skin conditions e.g.: psoriasis, eczema, dermatitis
- on open wounds, active infection or irritation of area of skin to be taped
- in case of cancer, kidney disease, congestive heart failure (CHF), deep venous thrombosis (DVT)
- if circulation compromised in area

6. WARNINGS & PRECAUTIONS

Keep out of reach of children. Do not re-use. Do not apply directly on wounds.

7. COMPLAINTS & CUSTOMER SERVICE SUPPORT

For questions, complaints or reporting an adverse reaction, you can contact the manufacturer.

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